


Customer information

Filtration · Rapid Tests · Water Analysis · Chromatography · Bioanalysis
Filtration · Schnellteste · Wasseranalytik · Chromatographie · Bioanalytik

Analytical Quality Control

July 2011

NANOCOLOR® NANOTURB



Method

Turbidity standard for nephelometric turbidity measurements for NANOCOLOR® spectrophotometers ($^{UV}_{VIS}$ und VIS)

Measuring range: 1 – 400 NTU

Shelf life: 8 months

REF: 925 702

List price: € 114.90 (valid until Dec. 31st 2011)

Content reagent set including:

1 tube with test solution 1 NTU
1 tube with test solution 4 NTU
1 tube with blank solution
1 instruction leaflet (DE, EN)

1 tube with test solution 100 NTU
1 tube with test solution 400 NTU
1 LOT certificate

Application

- Municipal waste water treatment plants
- Industrial waste water treatment plants
- Chemical industry

Potential clients

- All customers who are running nephelometric turbidity measurements (also when testing CSB)

Benefits

- Suitable as a primary standard for calibration and testing of the nephelometric turbidity unit according to ISO 7027
- NANOCOLOR® NANOTURB is free of hazardous or carcinogenic substances that have to be declared
- The test solutions come ready to use and just have to be inserted into the photometer. Dilution steps or contact with chemicals is avoided.

For more information on procedure and interferences, please consult the instructions. You can download the MSDS and instructions online from www.mn-net.com or request them directly from MACHEREY-NAGEL.

If you have further questions, please do not hesitate to contact us:

Technical support and customer service:

+49-(0)-2421 / 969-138
+49-(0)-2421 / 969-161
+49-(0)-2421 / 969-174
+49-(0)-2421 / 969-187
+49-(0)-2421 / 969-190

e-mail: csc@mn-net.com

To order:

e-mail: sales-de@mn-net.com

www.mn-net.com



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany
Germany and international:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

Switzerland:
MACHEREY-NAGEL AG
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

France:
MACHEREY-NAGEL EURL
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

USA:
MACHEREY-NAGEL Inc.
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com